



# **Plant Co-Immunoprecipitation (Co-IP) Kit**

**Catalog#JKR23001**

**Instruction Manual (For Two Groups)**

Sufficient reagents for 6T Co-IP assays per kit.


Store at -20 & 4°C

# Table of Contents

1.Experimental Principle .....	02
2.Technical Roadmap .....	02
3.Kit Components .....	04
4.Operating Steps .....	04
5.Frequently Asked Questions .....	06


## 1. Experimental Principle

Co-Immunoprecipitation (Co-IP) is a method based on the specific interaction between antigens and antibodies, used to study protein-protein interactions. After antibodies bind to their corresponding proteins in the lysate, they are incubated with Protein A/G-coupled Sepharose or Magnetic Beads. The Protein A/G bead-antibody-target protein complex is then obtained via centrifugation or a magnetic stand. Under high temperature and reducing agents, the antigen and antibody dissociate, and the supernatant—containing antibodies, target proteins, and a small amount of contaminating proteins—is collected. Proteins are subsequently identified by Western Blot or mass spectrometry (MS). Its schematic diagram is as follows:




## 2. Technical Roadmap

### 2.1 Experimental Flowchart



## 2.2 Control Setup Flowchart



### 3.Kit Components

Component	Volume (6T)	Storage Temperature
Plant lysis buffer	9mL	4°C
ProteinA/G Magnetic Beads	200 $\mu$ L	4°C
Incubation buffer	4.5mL	4°C
Wash buffer	20mL	4°C
Elution buffer	700 $\mu$ L	4°C, protect from light
Plant protease inhibitor (100 X)	35 $\mu$ L	-20°C
Normal Rabbit IgG (1mg/mL)	30 $\mu$ L	-20°C
Normal Mouse IgG (1mg/mL)	30 $\mu$ L	-20°C

Special reminder 1: Reagents and consumables such as PBS, 6X loading buffer, and magnetic rack need to be prepared by the user.

Special reminder 2: 6T refers to 6 single-group (1 IP group or 1 IgG group) immunoprecipitation experiments; the subsequent steps include 1 IgG group and 1 IP group, consuming 2T of reagent.

### 4.Operating Steps

#### 4.1 Total protein extraction

##### 4.1.1 Plant samples

- 1) Grinding: Take fresh young tissue or cryopreserved tissue (approximately 0.3g), place it in a sterilized and pre-cooled mortar, and grind with liquid nitrogen until powdered;
- 2) Lysis: Mix 990  $\mu$ L Plant Lysis buffer and 10  $\mu$ L Plant Protease inhibitor (100X) as the lysis solution, pipette 800  $\mu$ L of the lysis solution into the mortar, continue grinding on ice for 5-10 min until the sample forms a fine homogenate, transfer to a new EP tube, then add the remaining 200  $\mu$ L lysis solution to the mortar to collect residual sample, and similarly transfer it to the EP tube;
- 3) Fully lyse the EP tube containing the sample homogenate on ice for 30 min , inverting and mixing once every 5 min ;
- 4) Ultrasonication: Sonicate using an ultrasonic cell disruptor for 12–15 min at 20% power, with 3 sec of sonication and 3 sec of intervals, performed in an ice bath;
- 5) Centrifugation: 4°C, 12000 rpm, 10 min, collect the supernatant, then add Plant Lysis

buffer to the supernatant to bring the total volume to 1 mL, mix well.

Note: Perform the entire protein extraction process on ice to reduce protein degradation caused by high temperature; avoid bubble formation during sonication to minimize protein degradation. Store the total protein after lysis at -20°C.

## 4.2 CoIP

### 4.2.1 Antibody-bound Protein

- 1) Prepare two 1.5 mL EP tubes labeled as IP and IgG; add 3-5  $\mu\text{g}$  of the target antibody to the IP group, and add 3-5 $\mu\text{g}$  of the same species IgG to the IgG tube;
- 2) Take 100  $\mu\text{L}$  of the total protein solution from step 4.1, label it as Input, and store at -20°C for later use. Add 450  $\mu\text{L}$  of total protein solution to the IgG and IP tubes respectively, then supplement the volume to 800  $\mu\text{L}$  with Incubation buffer. Incubate with gentle mixing overnight (approximately 16 h) at 4°C.

### 4.2.2 Magnetic Bead Preparation

- 1) Remove the Protein A/G Magnetic Beads from the 4°C refrigerator, invert several times to mix the beads and solution thoroughly. Transfer 30 $\mu\text{L}$  each into two new 1.5 mL EP tubes, labeled as IgG and IP.
- 2) Add 0.5 mL of pre-chilled Wash buffer to both the IgG and IP tubes to resuspend the magnetic beads, place them on a magnetic stand for 1 min to separate the beads from the solution, carefully aspirate and discard the supernatant with a pipette, and repeat this step three times.

### 4.2.3 Magnetic Bead and Complex Binding

- 1) Add the two sets of incubated mixtures to the corresponding washed magnetic bead tubes, and incubate with gentle mixing at room temperature for 2 h;
- 2) Place the two tubes on a magnetic stand for 1 min to separate the beads from the solution, and carefully aspirate and discard the supernatant with a pipette.
- 3) Add 0.5 mL of pre-chilled Wash buffer to both the IgG and IP tubes, place them on a magnetic stand for 1min to separate the beads from the solution, carefully aspirate and discard the supernatant with a pipette, and repeat this step three times.

### 4.2.4 Elution

Add 100  $\mu\text{L}$  of Elution buffer to both the IgG and IP tubes, incubate in a boiling water bath for 10 min, place on a magnetic stand for 2 min, collect the supernatant as the eluate, mix with 20  $\mu\text{L}$  of 6 x SDS Loading buffer from the Input group, boil for 10 min, and store at -20°C for future use.

### 4.3 WB (Optional)

Take 20  $\mu$ L each of the IgG, IP, and Input group samples obtained in step 4.2.4 to perform WB detection.

### 4.4 Silver Staining (Optional)

- 1) Take 20  $\mu$ L each of the IgG, IP, and Input group samples obtained in step 4.2.4 to perform SDS-PAGE gel electrophoresis.
- 2) Fixation: Peel the post-electrophoresis gel from the glass plate, rinse it clean with water, place it in a clean 12 cm diameter glass dish, add deionized water to cover the gel, cover with a lid, shake at room temperature on a decolorizing shaker for 5 min, discard the deionized water, add fixation solution to cover the gel, cover with a lid, and shake at room temperature on a decolorizing shaker for 30 min.
- 3) Sensitization: Discard the fixative, add deionized water to cover the gel, cover with lid, shake at room temperature on a destaining shaker for 5 min, repeat washing with water once, for a total of 2 times. Add sensitization solution to cover the gel, cover with lid, shake at room temperature on a destaining shaker for 30 min.
- 4) Staining: Discard the sensitization solution, add deionized water to cover the gel, cover with lid, shake at room temperature on a destaining shaker for 2 min, repeat washing with water once, for a total of 2 times. Add staining solution to cover the gel, cover with lid, shake at room temperature on a destaining shaker for 20 min.
- 5) Color Development : Discard the staining solution, add deionized water to cover the gel, cover the lid, shake at room temperature on a destaining shaker for 1 min, repeat the water wash once, for a total of 2 times. Add developing solution to cover the gel, shake at room temperature on a destaining shaker for about 2 min until the solution turns cloudy. Discard the liquid, add fresh developing solution and continue developing until the target bands are clear, then photograph.

### 4.5 Mass Spectrometry (Optional)

Take 30  $\mu$ L of IgG and IP histone samples for LC-MS detection.

## 5. Frequently Asked Questions

### Q1: After CoIP followed by WB verification, no desired target band is found?

A: Caused by various reasons:

- 1) It is possible that the sample was degraded by proteases; the corresponding strategy is to add protease inhibitors, perform all operations on ice below 4°C, and avoid repeated freeze-thaw cycles.
- 2) It may be due to low antibody concentration resulting in faint bands; then adjust the IP or WB antibody concentration, and if necessary, set up a concentration gradient to determine

the optimal concentration.

- 3) The antibody affinity is too low; select an antibody suitable for IP or WB.
- 4) Some IP antibodies are not bound to the magnetic beads; in this case, choose magnetic beads suitable for IP.
- 5) If the Tag is not exposed on the surface of the fusion protein conformation, change the Tag fusion expression site.
- 6) The salt or alkalinity of the lysis buffer is too high; use a lysis buffer with lower salt or alkalinity. Inappropriate antibody selection; replace the antibody.

**Q2: Validation by WB after CoIP shows that although the target band is visible, the background is high:**

**A: Caused by multiple factors:**

- 1) High background due to non-specific protein binding. To avoid non-specific protein binding, lyse cells in serum- free solution, pre-wash immunoprecipitation with protein A/G magnetic beads, and increase rinse frequency and salt/alkalinity (high salt or detergent) after immunoprecipitation.
- 2) Contamination of experimental instruments or reagents; use clean instruments and reagents.
- 3) High background due to non-specific adsorption on the transfer membrane; wear gloves during experimental operations, use tweezers to handle, and avoid touching the membrane transfer surface.
- 4) If there may be large protein complexes that are not fully dissolved in the prepared sample, perform brief sonication after sample preparation (3 times, 5 sec each), then centrifuge and take the supernatant for subsequent experiments.
- 5) If washing is insufficient, perform multiple washes and consider increasing the concentration of NaCl and detergent in the washing solution.
- 6) If non-specific proteins may adsorb to the beads, perform preclearing to eliminate non-specific adsorption.
- 7) If the antibody itself has poor specificity, possibly leading to high background, select an appropriate antibody, considering monoclonal antibodies.
- 8) If using too many cells or tissues for lysis results in high background, reduce the sample amount, recommending 100-500  $\mu$ g of cell lysate.
- 9) Protein degradation may also result in high background; try to use freshly prepared samples whenever possible.



## **WUHAN GENECREATE BIOLOGICAL ENGINEERING CO.,LTD.**

**Address :** Wuhan Shengzhiyuan Biotechnology Innovation Industrial Park,  
North Shendun 5 th Road, Jiangxia District, Wuhan City, Hubei Province.

**Phone:** 027-87960366

**Email:** [marketing@genecreate.com](mailto:marketing@genecreate.com)

**Website:** [www.genecreate.cn](http://www.genecreate.cn)